

WARRANTY

The toothgrinder is covered by a 12 month warranty against defects in materials and fabrication (excluding battery). Within the warranty period, the toothgrinder will be repaired, adjusted or replaced. Return freight, insurance and proof of purchase must follow the product.

Damage caused by misuse is not covered by the warranty. There is no compensation for consequential damages. We exclude liability for indirect damage and loss.

WARNING

Remember that the cap should be in right place during grinding, groove is pressed around the tooth when grinding.

Uniteeth Toothgrinder Manual



INTRODUCTION

Uniteeth is a toothgrinder specially designed for grinding tusks of piglets. The grinding head consists of a fast rotating diamond coated cylinder and is constructed in a way so that the soft areas of the mouth is protected against damage.

Compared to tooth clipping with a tong, the toothgrinder is a better solution, both humanly and healthy, as the risks of diseases typically caused by tooth clipping by tong, is eliminated by using the toothgrinder instead.



Work smarter – not harder!

TOOTH GRINDING

Start the toothgrinder at the push button.

Place the grindstone on the tooth and press lightly. After approx 1 second 1/3 of the tooth has been removed. It is enough and best only to remove 1/3 of the tooth.

CHARGING

Remove the cable between the toothgrinder and the battery and connect to the charger instead.

Connect the charger with power supply. The red indicator lights up when battery is charging and turns off when the battery is fully charged.

The charger is a so-called two-stage charger, which means that the battery is not damaged by still charging even though it is fully charged. The charging function is suspended.

The battery is a lead storage battery type, by which the lifetime is reduced when it is not fully charged. As the charger at the same time do not damage the battery, it is best charge the battery after use, until the next time you have to use it.

- The battery has an effective worktime of approx. 45 minutes.
Similar to tooth grinding of approx. 100 piglets.
- The battery can with the included charger, not be overcharged.
- Charging time is approx. 10-14 hours.

CLEANING/DESINFECTION

The diamond head should be cleaned for tooth dust after use. Clean the diamond head with a damp toothbrush or similar. Disinfect the diamond head with rubbing alcohol or other disinfectant liquids.

REPLACING THE DIAMOND GRINDSTONE

By regular use, the diamond cylinder has a durability equal to grind approx. 7000 piglets. The diamond cylinder is changed by unscrewing the cap and loosen the cylinder with the provided 6 mm wrench.

MAINTENANCE

Uniteth should be stored in a dry environment, for example in an office or similar. Keep away from directly sunlight and temperatures over 40 degrees celcius. Remove regularly the diamond cylinder and clean the parts housed in the top. The parts should always be kept dry and clean to prevent rusting.

Charge the battery immediately after use.

TECHNICAL SPECIFICATIONS

Weight:	0.350 kg
Speed:	approx. 10.000 rpm
Charging time:	approx. 10-14 hours
Working time:	approx. 45 min.
Main supply:	230V
Volt:	12

